



APPLICATION NOTE

ICP - Mass Spectrometry

Authors:

Yongbo Dan

Honglan Shi

Xinhua Liang

Missouri University of
Science & Technology

Chady Stephan

PerkinElmer, Inc.
Shelton, CT

Measurement of Titanium Dioxide Nanoparticles in Sunscreen using Single Particle ICP-MS

Introduction

Titanium dioxide (TiO_2) nanoparticles are commonly used in sunscreens as part of

the formulation to block the sun's harmful UV rays. As the use of nanoparticles in consumer products has increased, concern has risen as to the health and environmental effects of nanoparticles. Since sunscreens contact skin and wash off in water, the TiO_2 nanoparticles can find their way into biological and environmental systems. As a result, there is a need to measure both the size and size distribution of these nanoparticles in sunscreens so as to assess their impact on human health and the environment.

Traditionally, nanoparticle size characteristics have been determined by several methods, including field flow fractionation (FFF), dynamic light scattering, and microscopy, among others. Recently, single particle inductively coupled plasma mass spectrometry (SP-ICP-MS) has been gaining attention as a way to both measure and characterize nanoparticles. Advantages of SP-ICP-MS include speed and the amount of information which can be gained from the analysis.

This work will focus on characterizing TiO₂ nanoparticles in sunscreens using SP-ICP-MS with a simple sample preparation.

Experimental

Sample Preparation

Sunscreen products were purchased at a local store and had the label information shown in Table 1. For each sample, 15 mL of sunscreen was added to a 50 mL container and mixed for three minutes to create a homogeneous sample. Next, 0.2 g of each homogenized sample was transferred to another container, followed by addition of 200 mL 1% Triton-X solution. These solutions were sonicated until all aggregates appeared to be broken up (5-10 minutes). After sonication, serial dilutions were performed on the samples with deionized water to produce particle concentrations of 100,000 – 200,000 particles/mL (as measured during analysis).

Transport efficiency determination was performed with gold nanoparticles (50, 80, 100 nm, nanoComposix, San Diego, California, U.S.A.) which were prepared by dilution in deionized water to a final nominal concentration of 100,000 particles/mL. All standards were sonicated for five minutes prior to analysis to ensure that any agglomerated particles were dispersed. Forty nanometer TiO₂ (US Research Nanomaterials Houston, Texas, U.S.A.) spikes were added to various samples to verify the accuracy of the TiO₂ size measurements.

Calibration standards for titanium consisted of 2, 5, and 10 ppb standards prepared in deionized water from a 1000 mg/L titanium standard (PerkinElmer, Shelton, Connecticut, U.S.A.).

Table 1. Label Information for Sunscreen Products.

| Sunscreen | SPF | TiO ₂ Content (%) |
|-----------|-----|------------------------------|
| 1 | 60+ | 4.9 |
| 2 | 50 | 6 |
| 3 | 45 | 0 |
| 4 | 50 | 6 |
| 5 | 45+ | 5.1 |

Instrument Conditions

All analyses were performed on a PerkinElmer NexION® 350D ICP-MS using the Syngistix™ Nano Application Software Module. Tables 2 and 3 show the instrumental and method parameters, respectively. When using the Syngistix Nano Application Module, the quadrupole settling time is automatically eliminated, thus ensuring that no particles are missed. The combination of elimination of quadrupole settling time and short dwell time results in the collection of multiple data points for each particle.

Table 2. NexION 350D ICP-MS Instrument Parameters.

| Parameter | Value |
|--------------------------|------------------|
| Nebulizer | Glass concentric |
| Spray chamber | Glass cyclonic |
| Sample uptake rate | 0.26-0.28 mL/min |
| RF power | 1600W |
| Analysis mode | Standard |
| Quadrupole settling time | 0 µs |

Table 3. Method Parameters.

| Parameter | Value |
|-----------------------|------------------------|
| Analyte | Ti48 |
| Dwell time | 100 µs |
| Data acquisition time | 100 sec |
| Density | 4.23 g/cm ³ |
| Ti mass fraction | 60% |

Results and Discussion

To characterize the Ti background of the system and reagents, a blank solution consisting of only the reagents (i.e. Triton X-100 in deionized water) was analyzed, with the resulting signal shown in Figure 1. While a few signal spikes are seen, the majority of the Ti background is below two counts, which demonstrates the cleanliness of the system.




Figure 1. Ti background signal for blank.

To determine the ability to see TiO₂ particles in a sunscreen matrix, the sunscreen without any TiO₂ particles was analyzed. Figure 2A shows the TiO₂ particle signal from a sunscreen which does not contain TiO₂ (according to its label). Aside from a few spikes, the background is equivalent to the blank (Figure 1).

Next, this sunscreen solution was spiked with 40 nm TiO₂ particles at a concentration of 6.65 µg/L; the resulting signal is shown in Figure 2B. Since each spike represents a particle, it is evident that TiO₂ particles can easily be seen in a sunscreen matrix.

Furthermore, the most frequent size particle was 37.7 µm, demonstrating the accuracy of the measurement.

Figure 3 shows the TiO₂ particle size distribution for a 20,000 times-diluted sunscreen (Sunscreen 5).

With the ability to accurately measure TiO₂ particles in sunscreen established, other sunscreen samples were analyzed three times; the results are shown in Table 4.





Table 4. TiO₂ Nanoparticle Analysis in Commercial Sunscreens.

| Sunscreen | Dilution Factor | Most Frequent Size (nm) | Particle Size Distribution (nm) | Particle Concentration (particles/mL) |
|-----------|-----------------|-------------------------|---------------------------------|---------------------------------------|
| 1 | 100,000 | 32 | 24 – 58 | 102,229 |
| 2 | 100,000 | 34 | 24 – 64 | 117,252 |
| 3* | 100,000 | - | - | - |
| 4 | 100,000 | 33 | 24 – 61 | 63,000 |
| 5 | 20,000 | 42 | 28 – 67 | 198,462 |

*This sunscreen contains no nanoparticles, as shown in Table 1.

These results indicate that the sunscreens had TiO₂ nanoparticles ranging from 24 – 67 nm, with similar size distributions in each. However, Sunscreen 5 is clearly different from the others: the smaller dilution factor required to obtain less than 200,000 particles/mL indicates that fewer TiO₂ particles are present in this sample than the others.

Conclusion

This work has demonstrated the ability to measure TiO₂ nanoparticles in commercial sunscreen samples. Using both the NexION 350 ICP-MS and the Syngistix Nano Application Software Module, the analysis is simple and rapid, yet can clearly differentiate the TiO₂ content among different samples.

Consumables Used

| Component | Description | Part Number |
|----------------------|--|-------------|
| Titanium standard | 1000 ppm, 125 mL | N9303806 |
| Sample uptake tubing | 0.38 mm id (green/orange), flared, 2-stop | N0777042 |
| Drain tubing | 1.30 mm id (gray/gray), Santoprene, 2-stop | N0777444 |