

Label Printing Capabilities:

LABWORKS LIMS and BarTender Software




www.perkinelmer.com

LABWORKS™ ES users can easily customize the software's paperwork settings to create the labels that best suit their laboratories' needs. One of the **LABWORKS™** ES label printing options is integration with BarTender® Label Printing Software from Seagull Scientific®.

As a **LABWORKS™** user, you can design labels as you see fit. The BarTender package gives you total control over type and size of text, bar codes and graphics. You can place objects exactly where needed, control rotation, import logos and images and export bar codes to other programs. In addition, individual text objects may contain multiple typefaces and styles.

The software also raises the power of rich text to a whole new level by allowing you to source different parts of a rich text object with external data fields.

All BarTender functionality can be controlled from within the **LABWORKS™** software. Departments from inventory and shipping to accounting and manufacturing can access BarTender's label printing capabilities. Within the **LABWORKS™** software, users can choose BarTender as the preferred label printing format. Users can then select data sources and specify network printers. **LABWORKS™** can also inquire about the current label status, including a Print Preview function so you can design quickly without wasting labels.



- ▶ Prompted Data Entry
- ▶ Predefined Industrial Symbols
- ▶ Microsoft® Office®-like User Interface
- ▶ Real-time On-screen Displays
- ▶ Industry-leading Data Access
- ▶ Over 20 International Languages
- ▶ Flexible Integration Across the Enterprise
- ▶ Supports any Database
- ▶ Multi-part Labels


FEATURES	BENEFITS
PROMPTED DATA ENTRY	Generate your own custom Data Entry screens to match your workflow.
PREDEFINED INDUSTRIAL SYMBOLS	Choose symbols from libraries such as electrical, hazardous material, packaging and more.
FAMILIAR MICROSOFT® OFFICE®-LIKE USER INTERFACE	Simple, intuitive mouse motions and easy-to-understand pop-up dialogs let you satisfy almost any compliance specification.
REAL-TIME ON-SCREEN DISPLAYS	Preview label design and layout in real-time without wasting labels.
INDUSTRY-LEADING DATA ACCESS	Read data from dozens of database formats, even from non-Windows platforms, and add native support for SAP IDoc files.
OVER 20 INTERNATIONAL LANGUAGES	Instantly switch BarTender's user interface into your preferred language.
FLEXIBLE INTEGRATION ACROSS THE ENTERPRISE	Allows non-Windows programs running on a variety of platforms to trigger and control label-printing functions. Advanced features include e-mailing of job status information across your enterprise and the Internet.
SUPPORTS ANY DATABASE	Simultaneously access multiple tables and data formats located on different computers.
MULTI-PART LABELS	Rotation of text and barcodes allows for different-sized label parts.

Security

Since it's often necessary to ensure that other users can't modify label designs, the software has a built-in security feature. The module offers Password Protection functionality by locking BarTender into a "print only" mode.

Flexible, Customizable Serialization

Generate anything from the simplest serial numbers to the most advanced, totally custom serialization sequences. You can specify numeric, alpha or a variety of predefined alphanumeric character sets or even specify completely customized sequences.



For more information, contact us at 800-762-4060 or visit us at www.labworks.com

PerkinElmer Life and Analytical Sciences
710 Bridgeport Avenue
Shelton, CT 06484-4794 USA
Phone: 800-762-4060 or
(+1) 203-925-4602
www.perkinelmer.com

