



DIRECT SAMPLE ANALYSIS

TOF MS


Food and Beverage: Caffeine in Decaffeinated Coffee

Easy to use

Product integrity

No chromatography

- Caffeine is a natural stimulant in tea and coffee.
- Decaffeinated coffee offers people the choice to avoid caffeine or reduce their caffeine intake
- DSA/TOF analysis of decaffeinated coffee confirmed caffeine was still present with high accurate mass
- Analysis was performed in 15 seconds confirmed product contamination with no sample preparation and external calibration



PerkinElmer, Inc.
940 Winter Street
Waltham, MA 02451 USA
P: (800) 762-4000 or
(+1) 203-925-4602
www.perkinelmer.com



For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs

Copyright © 2012, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.

010570_01