

Caution: For Laboratory Use. A product for research purposes only

**Sm-N1 labeled streptavidin**

**Product Number: AD0049**

Amount:	50 <i>ug</i>
Product lot no:	<input type="text" value="2621131"/>
Expiration date:	<input type="text" value="30 Aug 2020"/>
Labeling yield : (Eu / Protein)	<input type="text" value="7.78"/>
Concentration:	500 <i>ug/mL</i>
Buffer:	50 mM Tris-HCl buffered saline (pH 7.8) with azide as an antimicrobial agent and 0.1% BSA
Labeling reagent:	Sm-N1 chelate
Labeling conditions:	+4°C, pH 8.5, o/n
Purification method:	Gel filtration
Storage:	2-8 °C (36-46 °F)

For laboratory use. Reagents are stable, non-hazardous laboratory reagents and are not to be administered to humans or animals.

PerkinElmer, Inc.  
549 Albany Street  
Boston, MA 02118 USA  
P: (800) 762-4000 or (+1) 203-925-4602  
[www.perkinelmer.com/nenradiochemicals](http://www.perkinelmer.com/nenradiochemicals)  
**For a complete listing of our global offices, visit [www.perkinelmer.com/ContactUs](http://www.perkinelmer.com/ContactUs)**  
Copyright ©2010, PerkinElmer, Inc. All rights reserved.  
PerkinElmer® is a registered trademark of PerkinElmer, Inc.  
All other trademarks are the property of their respective owners.

