

Research Use Only. Not for use in diagnostic procedures.

Mouse IL10, lyophilized (0.3 µg)

Product No.: AL502S

Lot No.: 245-526-2

Material Provided

Format: 0.3 µg (one vial contains an amount of mouse IL10 sufficient for performing 10 AlphaLISA standard curves).
This product is designed to be used in combination with the AlphaLISA mIL10 Kit (AL502C or AL502F).

Manufacturing date: June 25, 2018

Product Information

Description: Mouse Interleukin 10 (mIL10), lyophilized (white powder).

Formulation: Lyophilized from a solution containing Tris, NaCl, Proclin-300, BSA, Dextran-10.

Storage: Store at +4°C.

Stability: This product is stable for at least 12 months from the manufacturing date when stored in its original packaging and the recommended storage conditions.

Reconstitution: Reconstitute mouse IL10 in 100 µL Milli-Q® grade H₂O and aliquot if required.

Aliquots should be stored in screw-capped polypropylene vials at -20°C. Reconstituted mouse IL10 is stable for at least 75 days at -20°C when stored in its original vial. Avoid multiple freeze-thaw cycles.

For more information regarding related AlphaLISA products and protocols refer to:

Technical Data Sheet: AlphaLISA mIL10 Kit, AL502C (500 assay points) & AL502F (5 000 assay points)

Technical Data Sheet: AlphaLISA Immunoassay Buffer (10X), AL000C (10 mL) & AL000F (100 mL)

Please visit our website www.perkinelmer.com/AlphaTech

This product is not for resale or distribution except by authorized distributors.

LIMITED WARRANTY: PerkinElmer BioSignal Inc. warrants that, at the time of shipment, the products sold by it are free from defects in material and workmanship and conform to specifications which accompany the product. PerkinElmer BioSignal Inc. makes no other warranty, express or implied with respect to the products, including any warranty of merchantability or fitness for any particular purpose. Notification of any breach of warranty must be made within 60 days of receipt unless otherwise provided in writing by PerkinElmer BioSignal Inc. No claim shall be honored if the customer fails to notify PerkinElmer BioSignal Inc. within the period specified. The sole and exclusive remedy of the customer for any liability of PerkinElmer BioSignal Inc. of any kind including liability based upon warranty (express or implied whether contained herein or elsewhere), strict liability contract or otherwise is limited to the replacement of the goods or the refunds of the invoice price of goods. PerkinElmer BioSignal Inc. shall not in any case be liable for special, incidental or consequential damages of any kind.

PerkinElmer, Inc.
940 Winter Street
Waltham, MA 02451 USA
P: (800) 762-4000 or
(+1) 203-925-4602
www.perkinelmer.com

