


Research Use Only. Not for use in diagnostic procedures.

**THYMIDINE, [METHYL-<sup>3</sup>H]-**

**Product Number: NET027Z**

**LOT SPECIFIC INFORMATION**

Lot Number:	202004
Specific Activity:	84 Ci/mmol
	3108 GBq/mmol
Production Date:	17-Mar-2020



M.W. 242  
C<sub>10</sub>H<sub>14</sub>N<sub>2</sub>O<sub>5</sub>

**PACKAGING:** 1 mCi/mL (37 MBq/ml) in aqueous solution, steri-packaged.

**STABILITY AND STORAGE RECOMMENDATIONS:**

When Thymidine, [methyl-<sup>3</sup>H]- is stored at refrigerated temperatures (5°C), in its original solvent and at its original concentration, the rate of decomposition is initially 4-5% per month. Stability is non-linear and not correlated to isotope half-life. Lot to lot variation may occur.

**SPECIFIC ACTIVITY RANGE:** 70-90 Ci/mmol (2590-3330 GBq/mmol)

**RADIOCHEMICAL PURITY:** This product was initially found to be greater than 97% when determined using the following method. The rate of decomposition can accelerate.

High pressure liquid chromatography on a C18 column using the following mobile phase:  
water : formic acid (conc.) : acetonitrile, (95:5:0.5).

**HAZARD INFORMATION:** WARNING: This product contains a chemical known to the state of California to cause cancer.

PerkinElmer, Inc.  
549 Albany Street  
Boston, MA 02118 USA  
P: (800) 762-4000 or (+1) 203-925-4602  
www.perkinelmer.com/enradiochemicals

For a complete listing of our global offices, visit [www.perkinelmer.com/ContactUs](http://www.perkinelmer.com/ContactUs)  
Copyright ©2010, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc.  
All other trademarks are the property of their respective owners.

