


Research Use Only. Not for use in diagnostic procedures.

THYMIDINE, [METHYL-³H]-

Product Number: NET027Z

LOT SPECIFIC INFORMATION

Lot Number:	202101
Specific Activity:	83.2 Ci/mmol
	3078 GBq/mmol
Production Date:	15-Dec-2020



M.W. 242
C₁₀H₁₄N₂O₅

PACKAGING: 1 mCi/mL (37 MBq/ml) in aqueous solution, steri-packaged.

STABILITY AND STORAGE RECOMMENDATIONS:

When Thymidine, [methyl-³H]- is stored at refrigerated temperatures (5°C), in its original solvent and at its original concentration, the rate of decomposition is initially 4-5% per month. Stability is non-linear and not correlated to isotope half-life. Lot to lot variation may occur.

SPECIFIC ACTIVITY RANGE: 70-90 Ci/mmol (2590-3330 GBq/mmol)

RADIOCHEMICAL PURITY: This product was initially found to be greater than 97% when determined using the following method. The rate of decomposition can accelerate.

High pressure liquid chromatography on a C18 column using the following mobile phase:
water : formic acid (conc.) : acetonitrile, (95:5:0.5).

HAZARD INFORMATION: WARNING: This product contains a chemical known to the state of California to cause cancer.

PerkinElmer, Inc.
549 Albany Street
Boston, MA 02118 USA
P: (800) 762-4000 or (+1) 203-925-4602
www.perkinelmer.com/enradiochemicals
For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs
Copyright ©2010, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc.
All other trademarks are the property of their respective owners.

